

[KO] CDK2 Mouse mAb

ABZ1066

Product Overview

Name	[KO] CDK2 Mouse mAb
Accession(Primary)	P24941
Host Species	Mouse
Isotype	IgG2a
Immunogen	Full length human recombinant protein of human CDK2 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB: 1/500-1/1000; ICC: 1/100-1/200; FC: 1/50-1/100

Target information (P24941)

Gene Name	CDK2
About Gene Symbol	CDK2 (Cyclin Dependent Kinase 2) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.