

[KO] CDK5 Mouse mAb

ABZ1067

Product Overview

Name	[KO] CDK5 Mouse mAb
Accession(Primary)	Q00535
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human CDK5 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB: 1/500-1/1000; ICC: 1/100-1/200

Target information (Q00535)

Gene Name	CDK5
About Gene Symbol	CDK5 (Cell division protein kinase 5) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.