

[KO] beta III Tubulin Mouse mAb

ABZ1069

Product Overview

Name	[KO] beta III Tubulin Mouse mAb
Accession(Primary)	Q13509
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human TUBB3 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 55 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (Q13509)

Gene Name	TUBB3
About Gene Symbol	TUBB3 (Tubulin Beta 3 Class III) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in fibrotic disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.