

[KO] YAP Mouse mAb

ABZ1070

Product Overview

Name	[KO] YAP Mouse mAb
Accession(Primary)	P46937
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human YAP1 protein fragments expressed in E.coli.
Calculated MW	Calculated MW: 54 kDa; Observed MW: 70 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (P46937)

Gene Name	YAP1
About Gene Symbol	YAP1 (Yes1 Associated Transcriptional Regulator) is a transcriptional coactivator involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.