

[KO] CDX2 Mouse mAb

ABZ1071

Product Overview

Name	[KO] CDX2 Mouse mAb
Accession(Primary)	Q99626
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human CDX2 produced in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q99626)

Gene Name	CDX2
About Gene Symbol	