

[KO] ZNRD2 Mouse mAb

ABZ1072

Product Overview

Name	[KO] ZNRD2 Mouse mAb
Accession(Primary)	O60232
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human ZNRD2 produced in 293T cell.
Calculated MW	Calculated MW: 21 kDa; Observed MW: 27 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; FC: 1/50-1/100

Target information (O60232)

Gene Name	ZNRD2
About Gene Symbol	