

[KO] p23 Mouse mAb

ABZ1073

Product Overview

Name	[KO] p23 Mouse mAb
Accession(Primary)	Q15185
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human PTGES3 produced in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 19 kDa; Observed MW: 23 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q15185)

Gene Name	PTGES3
About Gene Symbol	PTGES3 (Prostaglandin E Synthase 3) is a biosynthetic enzyme involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.