

[KO] GIPC1 Mouse mAb

ABZ1074

Product Overview

Name	[KO] GIPC1 Mouse mAb
Accession(Primary)	O14908
Host Species	Mouse
Isotype	IgG2a
Immunogen	Full length human recombinant protein of human GIPC1 produced in E.coli.
Calculated MW	Calculated MW: 36 kDa; Observed MW: 42 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (O14908)

Gene Name	GIPC1
About Gene Symbol	GIPC1 (GIPC PDZ Domain Containing Family Member 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechan