

[KO] SUGT1 Mouse mAb

ABZ1075

Product Overview

Name	[KO] SUGT1 Mouse mAb
Accession(Primary)	Q9Y2Z0
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human SUGT1 produced in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (Q9Y2Z0)

Gene Name	SUGT1
About Gene Symbol	