

[KO] HNRNPA0 Mouse mAb

ABZ1077

Product Overview

Name	[KO] HNRNPA0 Mouse mAb
Accession(Primary)	Q13151
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human HNRNPA0 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 31 kDa; Observed MW: 36 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q13151)

Gene Name	HNRNPA0
About Gene Symbol	HNRNPA0 (Heterogeneous Nuclear Ribonucleoprotein A0) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.