

## [KO] FERMT2 Mouse mAb

ABZ1078

### Product Overview

|                             |                                                                              |
|-----------------------------|------------------------------------------------------------------------------|
| <b>Name</b>                 | [KO] FERMT2 Mouse mAb                                                        |
| <b>Accession(Primary)</b>   | Q96AC1                                                                       |
| <b>Host Species</b>         | Mouse                                                                        |
| <b>Isotype</b>              | IgG2a                                                                        |
| <b>Immunogen</b>            | Full length human recombinant protein of human FERMT2 produced in 293T cell. |
| <b>Species Reactivity</b>   | Rat                                                                          |
| <b>Calculated MW</b>        | Calculated MW: 78 kDa; Observed MW: 78 kDa                                   |
| <b>Recommended Dilution</b> | WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/100-1/200                          |

### Target information (Q96AC1)

|                          |                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | FERMT2                                                                                                                                                                                                                                                          |
| <b>About Gene Symbol</b> | FERMT2 (FERM Domain Containing Kindlin 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to eluc |