

[KO] Chk2 Mouse mAb

ABZ1079

Product Overview

Name	[KO] Chk2 Mouse mAb
Accession(Primary)	O96017
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human CHEK2 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/100-1/200; FC: 1/50-1/100

Target information (O96017)

Gene Name	CHEK2
About Gene Symbol	CHEK2 (Checkpoint Kinase 2) is a serine/threonine protein kinase involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.