

[KO] CNPase Mouse mAb

ABZ1080

Product Overview

Name	[KO] CNPase Mouse mAb
Accession(Primary)	P09543
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human CNP produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P09543)

Gene Name	CNP
About Gene Symbol	CNP (2',3'-Cyclic Nucleotide 3' Phosphodiesterase) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.