

[KO] ADH5 Mouse mAb

ABZ1081

Product Overview

Name	[KO] ADH5 Mouse mAb
Accession(Primary)	P11766
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human ADH5 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P11766)

Gene Name	ADH5
About Gene Symbol	