

[KO] MYL6B Mouse mAb

ABZ1082

Product Overview

Name	[KO] MYL6B Mouse mAb
Accession(Primary)	P14649
Host Species	Mouse
Isotype	IgG2a
Immunogen	Full length human recombinant protein of human MYL6B produced in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (P14649)

Gene Name	MYL6B
About Gene Symbol	