

[KO] beta Catenin Mouse mAb

ABZ1083

Product Overview

Name	[KO] beta Catenin Mouse mAb
Accession(Primary)	P35222
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human CTNNB1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 100 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; FC: 1/50-1/100

Target information (P35222)

Gene Name	CTNNB1
About Gene Symbol	