

[KO] delta Catenin Mouse mAb

ABZ1084

Product Overview

Name	[KO] delta Catenin Mouse mAb
Accession(Primary)	O60716
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human CTNND1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 108 kDa; Observed MW: 108 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (O60716)

Gene Name	CTNND1
About Gene Symbol	CTNND1 (Catenin Delta 1) is a cell junction/adhesion protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.