

[KO] Lipoamide Dehydrogenase Mouse mAb

ABZ1085

Product Overview

Name	[KO] Lipoamide Dehydrogenase Mouse mAb
Accession(Primary)	P09622
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human DLD produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 40-70 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P09622)

Gene Name	DLD
About Gene Symbol	DLD (Dihydrolipoamide Dehydrogenase) is a metabolic dehydrogenase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and id