

[KO] Thymidine Phosphorylase Mouse mAb

ABZ1086

Product Overview

Name	[KO] Thymidine Phosphorylase Mouse mAb
Accession(Primary)	P19971
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human TYMP protein fragments expressed in E.coli.
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P19971)

Gene Name	TYMP
About Gene Symbol	