

[KO] EGFR Mouse mAb

ABZ1087

Product Overview

Name	[KO] EGFR Mouse mAb
Accession(Primary)	P00533
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human EGFR protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 134 kDa; Observed MW: 175 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P00533)

Gene Name	EGFR
About Gene Symbol	EGFR (Epidermal Growth Factor Receptor) is a receptor tyrosine kinase involved in immune signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.