

[KO] AK1 Mouse mAb

ABZ1088

Product Overview

Name	[KO] AK1 Mouse mAb
Accession(Primary)	P00568
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human AK1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (P00568)

Gene Name	AK1
About Gene Symbol	AK1 (Adenylate Kinase 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.