

[KO] AK3L1 Mouse mAb

ABZ1089

Product Overview

Name	[KO] AK3L1 Mouse mAb
Accession(Primary)	P27144
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human AK4 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 25 kDa; Observed MW: 30 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P27144)

Gene Name	AK4
About Gene Symbol	AK4 (Adenylate Kinase 4) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.