

[KO] ErbB 2 Mouse mAb

ABZ1090

Product Overview

Name	[KO] ErbB 2 Mouse mAb
Accession(Primary)	P04626
Host Species	Mouse
Isotype	IgG2b
Immunogen	Purified recombinant human ERBB2 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 138 kDa; Observed MW: 180 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P04626)

Gene Name	ERBB2
About Gene Symbol	ERBB2 (Erb-B2 Receptor Tyrosine Kinase 2) is a receptor tyrosine kinase involved in Hedgehog signaling and adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therape