

[KO] ERCC1 Mouse mAb

ABZ1091

Product Overview

Name	[KO] ERCC1 Mouse mAb
Accession(Primary)	P07992
Host Species	Mouse
Isotype	IgG2b
Immunogen	Full length human recombinant protein of human ERCC1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 40 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/100-1/200

Target information (P07992)

Gene Name	ERCC1
About Gene Symbol	ERCC1 (ERCC Excision Repair 1, Endonuclease Non-Catalytic Subunit) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.