

[KO] AKT1 Mouse mAb

ABZ1092

Product Overview

Name	[KO] AKT1 Mouse mAb
Accession(Primary)	P31749
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human AKT1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 56 kDa; Observed MW: 60 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; FC: 1/50-1/100

Target information (P31749)

Gene Name	AKT1
About Gene Symbol	AKT1 (AKT Serine/Threonine Kinase 1) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.