

[KO] ALDH3A2 Mouse mAb

ABZ1095

Product Overview

Name	[KO] ALDH3A2 Mouse mAb
Accession(Primary)	P51648
Host Species	Mouse
Isotype	IgG2b
Immunogen	Full length human recombinant protein of human ALDH3A2 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 50 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P51648)

Gene Name	ALDH3A2
About Gene Symbol	