

[KO] FGFR3 Mouse mAb

ABZ1096

Product Overview

Name	[KO] FGFR3 Mouse mAb
Accession(Primary)	P22607
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human FGFR3 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 88 kDa; Observed MW: 125 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P22607)

Gene Name	FGFR3
About Gene Symbol	FGFR3 (Fibroblast Growth Factor Receptor 3) is a receptor tyrosine kinase involved in adaptive immunity. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.