

[KO] FHL1 Mouse mAb

ABZ1097

Product Overview

Name	[KO] FHL1 Mouse mAb
Accession(Primary)	Q13642
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human FHL1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (Q13642)

Gene Name	FHL1
About Gene Symbol	