

[KO] FKBP5 Mouse mAb

ABZ1098

Product Overview

Name	[KO] FKBP5 Mouse mAb
Accession(Primary)	Q13451
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human FKBP5 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 51 kDa; Observed MW: 60 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q13451)

Gene Name	FKBP5
About Gene Symbol	