

[KO] TTLL12 Mouse mAb

ABZ1099

Product Overview

Name	[KO] TTLL12 Mouse mAb
Accession(Primary)	Q14166
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human TTLL12 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 74 kDa; Observed MW: 74 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/100-1/200; FC: 1/50-1/100

Target information (Q14166)

Gene Name	TTLL12
About Gene Symbol	