

[KO] SIRT5 Mouse mAb

ABZ1100

Product Overview

Name	[KO] SIRT5 Mouse mAb
Accession(Primary)	Q9NXA8
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human SIRT5 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q9NXA8)

Gene Name	SIRT5
About Gene Symbol	SIRT5 (Sirtuin 5) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets