

[KO] SIRT1 Mouse mAb

ABZ1101

Product Overview

Name	[KO] SIRT1 Mouse mAb
Accession(Primary)	Q96EB6
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human SIRT1 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 82 kDa; Observed MW: 110 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q96EB6)

Gene Name	SIRT1
About Gene Symbol	