

[KO] DDX58 Mouse mAb

ABZ1102

Product Overview

Name	[KO] DDX58 Mouse mAb
Accession(Primary)	O95786
Host Species	Mouse
Isotype	IgG2a
Immunogen	Full length human recombinant protein of human DDX58 produced in 293T cell.
Calculated MW	Calculated MW: 107 kDa; Observed MW: 107 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (O95786)

Gene Name	DDX58
About Gene Symbol	