

[KO] PITPNB Mouse mAb

ABZ1103

Product Overview

Name	[KO] PITPNB Mouse mAb
Accession(Primary)	P48739
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human PITPNB produced in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 32 kDa; Observed MW: 32 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P48739)

Gene Name	PITPNB
About Gene Symbol	