

[KO] PCSK9 Mouse mAb

ABZ1104

Product Overview

Name	[KO] PCSK9 Mouse mAb
Accession(Primary)	Q8NBP7
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human PCSK9 protein fragments expressed in E.coli.
Calculated MW	Calculated MW: 74 kDa; Observed MW: 65,80 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (Q8NBP7)

Gene Name	PCSK9
About Gene Symbol	