

[KO] SERBP1 Mouse mAb

ABZ1105

Product Overview

Name	[KO] SERBP1 Mouse mAb
Accession(Primary)	Q8NC51
Host Species	Mouse
Isotype	IgG2a
Immunogen	Purified recombinant human SERBP1 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 60 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q8NC51)

Gene Name	SERBP1
About Gene Symbol	