

[KO] TIPRL Mouse mAb

ABZ1106

Product Overview

Name	[KO] TIPRL Mouse mAb
Accession(Primary)	O75663
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human TIPRL produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 31 kDa; Observed MW: 35 kDa
Recommended Dilution	WB: 1/500-1/1000; FC: 1/50-1/100

Target information (O75663)

Gene Name	TIPRL
About Gene Symbol	