

[KO] GBA Mouse mAb

ABZ1107

Product Overview

Name	[KO] GBA Mouse mAb
Accession(Primary)	P04062
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human GBA protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60-70 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P04062)

Gene Name	GBA
About Gene Symbol	