

[KO] NPTN Mouse mAb

ABZ1108

Product Overview

Name	[KO] NPTN Mouse mAb
Accession(Primary)	Q9Y639
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human NPTN produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 55 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q9Y639)

Gene Name	NPTN
About Gene Symbol	