

[KO] TRMT2A Mouse mAb

ABZ1109

Product Overview

Name	[KO] TRMT2A Mouse mAb
Accession(Primary)	Q8IZ69
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human TRMT2A produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 69 kDa; Observed MW: 69 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q8IZ69)

Gene Name	TRMT2A
About Gene Symbol	TRMT2A (TRNA Methyltransferase 2 Homolog) is a methyltransferase involved in translational regulation. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms an