

[KO] FOXP1 Mouse mAb

ABZ1110

Product Overview

Name	[KO] FOXP1 Mouse mAb
Accession(Primary)	Q9H334
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human FOXP1 produced in 293T cell.
Calculated MW	Calculated MW: 75 kDa; Observed MW: 75 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q9H334)

Gene Name	FOXP1
About Gene Symbol	