

[KO] AMF Mouse mAb

ABZ1111

Product Overview

Name	[KO] AMF Mouse mAb
Accession(Primary)	P06744
Host Species	Mouse
Isotype	IgG2b
Immunogen	Full length human recombinant protein of human GPI produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 63 kDa; Observed MW: 63 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; FC: 1/50-1/100

Target information (P06744)

Gene Name	GPI
About Gene Symbol	GPI (Glucose-6-Phosphate Isomerase) is a metabolic enzyme involved in adaptive immunity. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.