

[KO] Granulin Mouse mAb

ABZ1112

Product Overview

Name	[KO] Granulin Mouse mAb
Accession(Primary)	P28799
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human GRN protein fragments expressed in E.coli.
Calculated MW	Calculated MW: 64 kDa; Observed MW: 70 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P28799)

Gene Name	GRN
About Gene Symbol	GRN (Granulin Precursor) is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.