

[KO] Glucocorticoid Receptor Mouse mAb

ABZ1113

Product Overview

Name	[KO] Glucocorticoid Receptor Mouse mAb
Accession(Primary)	P04150
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human NR3C1 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 86 kDa; Observed MW: 90-120 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P04150)

Gene Name	NR3C1
About Gene Symbol	NR3C1 (Nuclear Receptor Subfamily 3 Group C Member 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms.