

## [KO] PACSIN3 Mouse mAb

ABZ1115

### Product Overview

|                             |                                                                               |
|-----------------------------|-------------------------------------------------------------------------------|
| <b>Name</b>                 | [KO] PACSIN3 Mouse mAb                                                        |
| <b>Accession(Primary)</b>   | Q9UKS6                                                                        |
| <b>Host Species</b>         | Mouse                                                                         |
| <b>Isotype</b>              | IgG1                                                                          |
| <b>Immunogen</b>            | Full length human recombinant protein of human PACSIN3 produced in 293T cell. |
| <b>Species Reactivity</b>   | Rat                                                                           |
| <b>Calculated MW</b>        | Calculated MW: 48 kDa; Observed MW: 48 kDa                                    |
| <b>Recommended Dilution</b> | WB: 1/500-1/1000; IHC: 1/50-1/100                                             |

### Target information (Q9UKS6)

|                          |                                                                                                                                                                                                                                                                                                                                                      |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | PACSIN3                                                                                                                                                                                                                                                                                                                                              |
| <b>About Gene Symbol</b> | PACSIN3 (Protein Kinase C And Casein Kinase Substrate In Neurons 3) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |