

[KO] CERS2 Mouse mAb

ABZ1116

Product Overview

Name	[KO] CERS2 Mouse mAb
Accession(Primary)	Q96G23
Host Species	Mouse
Isotype	IgG2a
Immunogen	Purified recombinant human CERS2 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 38 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q96G23)

Gene Name	CERS2
About Gene Symbol	