

[KO] hnRNP F Mouse mAb

ABZ1119

Product Overview

Name	[KO] hnRNP F Mouse mAb
Accession(Primary)	P52597
Host Species	Mouse
Isotype	IgG2a
Immunogen	Full length human recombinant protein of human HNRNPF produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/100-1/200; FC: 1/50-1/100

Target information (P52597)

Gene Name	HNRNPF
About Gene Symbol	HNRNPF (Heterogeneous Nuclear Ribonucleoprotein F) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.