

[KO] Hsp40 Mouse mAb

ABZ1121

Product Overview

Name	[KO] Hsp40 Mouse mAb
Accession(Primary)	P25685
Host Species	Mouse
Isotype	IgG2a
Immunogen	Full length human recombinant protein of human DNAJB1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; FC: 1/50-1/100

Target information (P25685)

Gene Name	DNAJB1
About Gene Symbol	DNAJB1 (DnaJ Heat Shock Protein Family (Hsp40) Member B1) is a molecular chaperone involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.