

[KO] ING2 Mouse mAb

ABZ1124

Product Overview

Name	[KO] ING2 Mouse mAb
Accession(Primary)	Q9H160
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human ING2 produced in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 37 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (Q9H160)

Gene Name	ING2
About Gene Symbol	