

[KO] IRF3 Mouse mAb

ABZ1125

Product Overview

Name	[KO] IRF3 Mouse mAb
Accession(Primary)	Q14653
Host Species	Mouse
Isotype	IgG2b
Immunogen	Full length human recombinant protein of human IRF3 produced in 293T cell.
Calculated MW	Calculated MW: 47 kDa; Observed MW: 56 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q14653)

Gene Name	IRF3
About Gene Symbol	IRF3 (Interferon Regulatory Factor 3) is an interferon regulatory transcription factor involved in JAK-STAT/cytokine signaling and innate immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.